**ANCA - ETHANOL**

Indirect Immunofluorescence HUMAN NEUTROPHILS

**A.MENARINI**  
diagnostics

COD 50053 10 x 6 tests

COD 50054 10 x 12 tests

Only for *in vitro* use in the clinical laboratory**INTENDED USE**

Reagents for the qualitative determination of anti-neutrophil cytoplasmic antibodies (ANCA) in human serum by indirect immunofluorescence.

**PRINCIPLE OF THE METHOD**

Serum anti-neutrophil cytoplasmic antibodies (ANCA) bind to the corresponding antigens present in human neutrophils. The resulting antigen-antibody complexes are detected by means of a fluorescein labeled anti-human immunoglobulin, and visualized with the aid of a fluorescence microscope<sup>1</sup>.

**CONTENTS**

|                  |                         | COD 50053    | COD 50054     |
|------------------|-------------------------|--------------|---------------|
| <b>SLIDE</b>     | Slides                  | 10 x 6 tests | 10 x 12 tests |
| <b>PBS</b>       | PBS (10x)               | 1 x 100 mL   | -             |
| <b>CONTROL +</b> | P-ANCA Positive Control | 1 x 1 mL     | -             |
| <b>CONTROL +</b> | C-ANCA Positive Control | 1 x 1 mL     | -             |
| <b>CONTROL -</b> | Negative Control        | 1 x 1 mL     | -             |
| <b>CONJUGATE</b> | IgG FITC/Evans          | 1 x 5 mL     | -             |
| <b>MOUNT</b>     | Mounting Medium         | 1 x 5 mL     | -             |
| <b>COVER SLD</b> | Coverslips              | 1 x 12       | -             |

**COMPOSITION**

**SLIDE** Ethanol-fixed human neutrophils immobilized on each well.

**PBS** Sodium phosphate 112.5 mmol/L, potassium phosphate 30 mmol/L, sodium chloride 1.15 mol/L, sodium azide 0.95 g/L, pH 7.2.

**CONTROL+** **P-ANCA:** Human serum containing anti-neutrophil cytoplasmic antibodies (ANCA) "perinuclear pattern", sodium azide 0.95 g/L.

**CONTROL+** **C-ANCA:** Human serum containing anti-neutrophil cytoplasmic antibodies (ANCA) "cytoplasmic pattern", sodium azide 0.95 g/L.

**CONTROL -** Human serum, sodium azide 0.95 g/L.

**CONJUGATE** Goat anti-human IgG conjugated with fluorescein isothiocyanate (FITC), Evans blue 0.01 g/L, sodium azide 0.95 g/L.

**MOUNT** Mowiol 1%, Glycerol 87%, Tris 20 mmol/L, sodium azide 0.95 g/L.

*Human sera used in the preparation of the positive and negative controls have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, the controls should be handled cautiously as potentially infectious.*

**STORAGE**

Store at 2-8°C.

Reagents are stable until the expiry date shown on the label when stored tightly closed and if contamination is prevented during their use.

**Indications of deterioration:**

- Liquid components: Presence of particulate material, turbidity.
- Slides: rips in the sealing bag, macroscopic defects on the cell culture like scratches or monolayer peeling off.

**AUXILIARY REAGENTS**

**PBS** PBS (10x)

**CONJUGATE** IgG FITC/Evans

**MOUNT** Mounting Medium

**CONTROL + P-ANCA:** Human serum containing anti-neutrophil cytoplasmic antibodies (ANCA) "cytoplasmic pattern", sodium azide 0.95 g/L.

**CONTROL + C-ANCA:** Human serum containing antinuclear antibodies (ANA) speckled pattern, sodium azide 0.95 g/L.

**CONTROL –** Negative Control.

**REAGENT PREPARATION**

**PBS:** Dilute PBS 1/10 with distilled water. Stable 1 month at 2-8°C, if stored at the recommended temperature, well closed and care is taken to prevent contamination during its use.

All other reagents are provided ready to use.

**ADDITIONAL EQUIPMENT**

- Moist chamber
- Wash tray
- Fluorescence microscope equipped with a 495 nm excitation filter and a 525 nm emission filter for FITC visualization

**SAMPLES**

Serum or plasma collected by standard procedures. Stable 7 days at 2-8°C.

Dilute samples 1/20 in PBS (see Reagent Preparation) before the assay.

For titration of positive samples, make two-fold serial dilutions starting from 1/40 in PBS.

However, each laboratory should establish its own sample dilution based on the population characteristics.


**PROCEDURE**

1. Bring the reagents and samples to room temperature.
2. Place 1 drop (25  $\mu$ L) of the diluted sample or Control on each slide well, making sure that it is completely covered (Note 1).
3. Incubate the slide for 30 minutes at room temperature (15-30°C) into a moist chamber.
4. Drain sample drops off by gently tapping the inclined slide. Avoid cross-contamination of the sera.
5. Rinse gently the slide with PBS (see Reagent Preparation). (Note 2).
6. Wash thoroughly the slide by immersing in a washing tray filled with PBS for 5 minutes. Change PBS and repeat wash.
7. Carefully dry off the slides by using the blotting paper provided. Keep the substrate moist along the procedure.
8. Place 1 drop of Conjugate on each well. Incubate the slide for 30 minutes at room temperature (15-30°C) into a moist chamber.
9. Wash (step 6) and dry (step 7).

10. Place several drops of Mounting Medium on the slide and cover with a coverslip avoiding the formation of air bubbles.

## READING

Examine the slide using the fluorescence microscope (250-400x). For best results, the slides should be read immediately. Select reading fields in the area indicated in the picture, between the center area and the edge area. Select fields with uniform spacing between cells and uniform nuclei brightness. Fluorescent intensity in the center and edge areas is not representative of the slide preparation.

Observation of specific fluorescent staining, described as follows at the recommended dilution, should be considered as a positive result.

There are two main staining patterns in **ethanol**-fixed neutrophils associated with autoimmune vasculitis: a diffuse granular cytoplasmic staining with higher interlobular intensity (C-ANCA), and a compact staining of the perinuclear zone of the cytoplasm (P-ANCA).

In some cases, the P-ANCA staining pattern observed on ethanol-fixed neutrophils may be due to antinuclear antibodies (ANA). This result should be confirmed by running the sample on **formalin**-fixed neutrophils.

There is only one staining pattern on formalin-fixed neutrophils associated with autoimmune vasculitis: C-ANCA

In addition to these staining patterns, other less frequent patterns, like X-ANCA (or atypical P-ANCA), can be identified, showing a broad non homogeneous perinuclear staining, eventually accompanied by a diffuse cytoplasmic pattern with no accentuation of the interlobular zone. For a proper observation of this pattern, the use of **methanol**-fixed neutrophils is recommended.

Positive sera may be tittered. The serum titer is defined as the highest dilution showing a positive result.

When none of the above specific stainings are observed, the result should be considered negative for these autoantibodies.

## QUALITY CONTROL

Positive Controls (CONTROL +) and Negative Control (CONTROL -) provided with kits cod 50053 should be tested together with the patients samples, in order to verify the assay performance.

Positive Controls (CONTROL +) should give the above described specific staining.

Negative Control (CONTROL -) should not give any specific staining.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

## ASSAY CHARACTERISTICS

The IgG FITC/Evans conjugate is calibrated versus the International Standard of the WHO for sheep anti-human immunoglobulins conjugated with FITC.

The specificity of the C-ANCA Positive Control has been verified against the PR3-ANCA Human Reference serum #16 from the CDC. The specificity of the P-ANCA Positive Control has been verified against the PR3-ANCA Human Reference serum #15 from the CDC.

The C-ANCA staining pattern on ethanol-fixed neutrophils is the result of the antibody binding mainly to the antigen proteinase 3. In formalin-fixed neutrophils, this staining pattern is due to the antibody binding to myeloperoxidase.

The P-ANCA staining pattern on ethanol-fixed neutrophils is the result of the antibody binding to myeloperoxidase.

### **DIAGNOSTIC CHARACTERISTICS**

There is a strong association between Wegener's granulomatosis and the presence of C-ANCA pattern associated with PR-3. This staining pattern is almost exclusively found in patients with this disease and its diagnostic specificity is higher than 90%. Sensitivity depends on the disease stage but can be found in average around 66%<sup>(2)</sup>. P-ANCA pattern associated with MPO can be found in necrotizing idiopathic vasculitis, like microscopic polyangiitis, idiopathic crescentic glomerulonephritis, Churg-Strauss syndrome, polyarteritis nodosa and Wegener's granulomatosis. Sensitivities of P-ANCA for microscopic polyangiitis and Churg-Strauss syndrome are 58% and 74.5%, respectively, whereas specificities are 81% and 95%<sup>(3,4)</sup>. X-ANCA pattern is associated with the diagnostic of inflammatory bowel diseases like ulcerative colitis (50-70%) and Crohn's disease (10-30%), as well as other nonvasculitic autoimmune diseases<sup>(5)</sup>.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

### **NOTES**

1. Avoid touching the cells fixed into the wells along the procedure.
2. Use a squeeze bottle or a pipette to wash the slides, avoiding cross-contamination among the adjacent samples.

### **BIBLIOGRAPHY**

1. Melnickoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Rao JK, Weinberger M, Oddone EZ, Allen NB, Landsman P, Feussner JR The Role of Antineutrophil Cytoplasmic Antibody (c-ANCA) Testing in the Diagnosis of Wegener Granulomatosis: A Literature Review and Meta-analysis. *Ann Intern Med.* 1995;123(12):925-932.
3. Hagen EC, Daha MR, Hermans J, Andrassy K, Csernok E, Gaskin G, Lesavre P, Lüdemann J, Rasmussen N, Sinico RA, Wiik A, van der Woude FJ. Diagnostic value of standardized assays for anti-neutrophil cytoplasmic antibodies in idiopathic systemic vasculitis. EC/BCR Project for ANCA Assay Standardization. *Kidney Int.* 1998; 53 (3): 743-53
4. Sinico RA, Di Toma L, Maggiore U, Bottero P, Radice A, Tosoni C, Grasselli C, Pavone L, Gregorini G, Monti S, Frassi M, Vecchio F, Corace C, Venegoni E, Buzio C. Prevalence and clinical significance of antineutrophil cytoplasmic antibodies in Churg-Strauss syndrome. *Arthritis Rheum.* 2005;52(9):2589-93.
5. Savage J, Dimech W, Fritzier M, et al. Addendum to the international Consensus Statement on Testing and Reporting of Antineutrophil Cytoplasmic Antibodies. *Am J Clin Pathol* 2003;120:312-318.

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## **BIBLIOGRAFIA**

1. Melnickoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Rao JK, Weinberger M, Oddone EZ, Allen NB, Landsman P, Feussner JR The Role of Antineutrophil Cytoplasmic Antibody (c-ANCA) Testing in the Diagnosis of Wegener Granulomatosis: A Literature Review and Meta-analysis. Ann Intern Med. 1995;123(12):925-932.
3. Hagen EC, Daha MR, Hermans J, Andrassy K, Csernok E, Gaskin G, Lesavre P, Lüdemann J, Rasmussen N, Sinico RA, Wiik A, van der Woude FJ. Diagnostic value of standardized assays for antineutrophil cytoplasmic antibodies in idiopathic systemic vasculitis. EC/BCR Project for ANCA Assay Standardization. Kidney Int. 1998; 53 (3): 743-53
4. Sinico RA, Di Toma L, Maggiore U, Bottero P, Radice A, Tosoni C, Grasselli C, Pavone L, Gregorini G, Monti S, Frassi M, Vecchio F, Corace C, Venegoni E, Buzio C. Prevalence and clinical significance of antineutrophil cytoplasmic antibodies in Churg-Strauss syndrome. Arthritis Rheum. 2005;52(9):2589-93.
5. Savige J, Dimech W, Fritzler M, et al. Addendum to the international Consensus Statement on Testing and Reporting of Antineutrophil Cytoplasmic Antibodies. Am J Clin Pathol 2003;120:312-318.

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1. Melnickoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Rao JK, Weinberger M, Oddone EZ, Allen NB, Landsman P, Feussner JR The Role of Antineutrophil Cytoplasmic Antibody (c-ANCA) Testing in the Diagnosis of Wegener Granulomatosis: A Literature Review and Meta-analysis. *Ann Intern Med.* 1995;123(12):925-932.
3. Hagen EC, Daha MR, Hermans J, Andrassy K, Csernok E, Gaskin G, Lesavre P, Lüdemann J, Rasmussen N, Sinico RA, Wiik A, van der Woude FJ. Diagnostic value of standardized assays for antineutrophil cytoplasmic antibodies in idiopathic systemic vasculitis. EC/BCR Project for ANCA Assay Standardization. *Kidney Int.* 1998; 53 (3): 743-53
4. Sinico RA, Di Toma L, Maggiore U, Bottero P, Radice A, Tosoni C, Grasselli C, Pavone L, Gregorini G, Monti S, Frassi M, Vecchio F, Corace C, Venegoni E, Buzio C. Prevalence and clinical significance of antineutrophil cytoplasmic antibodies in Churg-Strauss syndrome. *Arthritis Rheum.* 2005;52(9):2589-93.
5. Savige J, Dimech W, Fritzler M, et al. Addendum to the international Consensus Statement on Testing and Reporting of Antineutrophil Cytoplasmic Antibodies. *Am J Clin Pathol* 2003;120:312-318.

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## BIBLIOGRAPHIE

1. Melnicoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Rao JK, Weinberger M, Oddone EZ, Allen NB, Landsman P, Feussner JR. The Role of Antineutrophil Cytoplasmic Antibody (c-ANCA) Testing in the Diagnosis of Wegener Granulomatosis: A Literature Review and Meta-analysis. *Ann Intern Med.* 1995;123(12):925-932.
3. Hagen EC, Daha MR, Hermans J, Andrassy K, Csernok E, Gaskin G, Lesavre P, Lüdemann J, Rasmussen N, Sinico RA, Wiik A, van der Woude FJ. Diagnostic value of standardized assays for antineutrophil cytoplasmic antibodies in idiopathic systemic vasculitis. EC/BCR Project for ANCA Assay Standardization. *Kidney Int.* 1998; 53 (3): 743-53
4. Sinico RA, Di Toma L, Maggiore U, Bottero P, Radice A, Tosoni C, Grasselli C, Pavone L, Gregorini G, Monti S, Frassi M, Vecchio F, Corace C, Venegoni E, Buzio C. Prevalence and clinical significance of antineutrophil cytoplasmic antibodies in Churg-Strauss syndrome. *Arthritis Rheum.* 2005;52(9):2589-93.
5. Savige J, Dimech W, Fritzer M, et al. Addendum to the international Consensus Statement on Testing and Reporting of Antineutrophil Cytoplasmic Antibodies. *Am J Clin Pathol* 2003;120:312-318.

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1. Melnickoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Rao JK, Weinberger M, Oddone EZ, Allen NB, Landsman P, Feussner JR The Role of Antineutrophil Cytoplasmic Antibody (c-ANCA) Testing in the Diagnosis of Wegener Granulomatosis: A Literature Review and Meta-analysis. Ann Intern Med. 1995;123(12):925-932.
3. Hagen EC, Daha MR, Hermans J, Andrassy K, Csernok E, Gaskin G, Lesavre P, Lüdemann J, Rasmussen N, Sinico RA, Wiik A, van der Woude FJ. Diagnostic value of standardized assays for antineutrophil cytoplasmic antibodies in idiopathic systemic vasculitis. EC/BCR Project for ANCA Assay Standardization. Kidney Int. 1998; 53 (3): 743-53
4. Sinico RA, Di Toma L, Maggiore U, Bottero P, Radice A, Tosoni C, Grasselli C, Pavone L, Gregorini G, Monti S, Frassi M, Vecchio F, Corace C, Venegoni E, Buzio C. Prevalence and clinical significance of antineutrophil cytoplasmic antibodies in Churg-Strauss syndrome. Arthritis Rheum. 2005;52(9):2589-93.
5. Savige J, Dimech W, Fritzer M, et al. Addendum to the international Consensus Statement on Testing and Reporting of Antineutrophil Cytoplasmic Antibodies. Am J Clin Pathol 2003;120:312-318.

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## LITERATUR

1. Melnicoff MJ. Immunfluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Rao JK, Weinberger M, Oddone EZ, Allen NB, Landsman P, Feussner JR The Role of Antineutrophil Cytoplasmic Antibody (c-ANCA) Testing in the Diagnosis of Wegener Granulomatosis: A Literature Review and Meta-analysis. Ann Intern Med. 1995;123(12):925-932.
3. Hagen EC, Daha MR, Hermans J, Andrassy K, Csernok E, Gaskin G, Lesavre P, Lüdemann J, Rasmussen N, Sinico RA, Wiik A, van der Woude FJ. Diagnostic value of standardized assays for antineutrophil cytoplasmic antibodies in idiopathic systemic vasculitis. EC/BCR Project for ANCA Assay Standardization. Kidney Int. 1998; 53 (3): 743-53
4. Sinico RA, Di Toma L, Maggiore U, Bottero P, Radice A, Tosoni C, Grasselli C, Pavone L, Gregorini G, Monti S, Frassi M, Vecchio F, Corace C, Venegoni E, Buzio C. Prevalence and clinical significance of antineutrophil cytoplasmic antibodies in Churg-Strauss syndrome. Arthritis Rheum. 2005;52(9):2589-93.
5. Savige J, Dimech W, Fritzler M, et al. Addendum to the international Consensus Statement on Testing and Reporting of Antineutrophil Cytoplasmic Antibodies. Am J Clin Pathol 2003;120:312-318.

**MANUFACTURER:** BioSystems S.A.  
**PRODUTTORE:** Costa Brava 30  
**FABRICANTE:** 08030 Barcelona (Spain)  
**FABRICANT:** Tel: +34 93 311 00 00  
**KATAKEYAZTHI:** Fax: +34 93 311 76 09  
**HERSTELLER:** Mail: Biosystems@biosystems.es

**DISTRIBUTOR:** A. Menarini Diagnostics SRL  
**DISTRIBUTORE:** Via Sette Santi 3  
**DISTRIBUIDOR:** 50131 Firenze  
**DISTRIBUTEUR:** Tel. 0039-055-56801  
**ΔΙΑΝΟΜΕΑΣ:** Fax 0039-055-5680-902  
**HÄNDLER:**