AUTOANTIBODIES rKSL

Indirect Immunofluorescence RAT KIDNEY/STOMACH/LIVER

A.MENARINI
diagnostics

COD 50065 12 x 8 tests

COD 50052 12 x 8 tests

Only for *in vitro* use in the clinical laboratory

INTENDED USE

Reagents for the qualitative determination of autoimmune antibodies in human serum by indirect immunofluorescence. The obtained results are useful as an aid in the diagnosis of several organ specific autoimmune diseases.

PRINCIPLE OF THE METHOD

Serum autoantibodies directed against the nucleus (ANA), the mitochondria (AMA), the smooth muscle (ASMA), the gastric parietal cell (APCA), the liver-kidney microsomes (LKM), the reticulin and others, bind to the corresponding antigens present in a preparation containing rat liver, kidney and stomach. The resulting antigen-antibody complexes are detected by means of a fluorescein labeled anti-human immunoglobulin, and visualized with the aid of a fluorescence microscope¹.

CONTENTS

		COD 50065	COD 50052
SLIDE	Slides	12 x 8 tests	12 x 8 tests
PBS	PBS (10x)	1 x 100 mL	-
CONTROL +	AMA Positive Control	1 x 1 mL	-
CONTROL -	Negative Control	1 x 1 mL	-
CONJUGATE	IgG FITC/Evans	1 x 5 mL	-
MOUNT	Mounting Medium	1 x 5 mL	-
COVER SLD	Coverslips	1 x 12	-

COMPOSITION

SLIDE Sections of rat liver, kidney and stomach.

PBS Sodium phosphate 112.5 mmol/L, potassium phosphate 30 mmol/L, sodium chloride 1.15 mol/L, sodium azide 0.95 g/L, pH 7.2.

CONTROL+ Human serum containing anti-mitochondrial antibodies (AMA), sodium azide 0.95 g/L.

CONTROL - Human serum, sodium azide 0.95 g/L.

CONJUGATE Goat anti-human IgG conjugated with fluorescein isothiocyanate (FITC), Evans blue 0.01 g/L, sodium azide 0.95 g/L.

MOUNT Mowiol 1%, Glycerol 87%, Tris 20 mmol/L, sodium azide 0.95 g/L.

Human sera used in the preparation of the positive and negative controls have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, the controls should be handled cautiously as potentially infectious.

STORAGE

Store at 2-8°C.

Reagents are stable until the expiry date shown on the label when stored tightly closed and if contamination is prevented during their use.

Indications of deterioration:

- Liquid components: Presence of particulate material, turbidity.
- Slides: rips in the sealing bag, macroscopic defects on the cell culture like scratches or monolayer peeling off.

AUXILIARY REAGENTS

PBS PBS (10x).

CONJUGATE IgG FITC/Evans.

MOUNT Mounting Medium.

CONTROL + ANA-Ho Positive Control: Human serum containing anti-nuclear antibodies (ANA) homogeneous pattern, sodium azide 0.95 g/L.

CONTROL + AMA Positive Control: Human serum containing anti-mitochondrial antibodies (AMA), sodium azide 0.95 g/L.

CONTROL + ASMA Positive Control: Human serum containing anti-smooth muscle antibodies (ASMA), sodium azide 0.95 g/L.

CONTROL + LKM Positive Control: Human serum containing anti-liver-kidney microsomal antibodies (LKM), sodium azide 0.95 g/L.

CONTROL + APCA Positive Control: Human serum containing anti-parietal cell muscle antibodies (APCA), sodium azide 0.95 g/L.

CONTROL – Negative Control.

REAGENT PREPARATION

PBS: Dilute PBS 1/10 with distilled water. Stable 1 month at 2-8°C, if stored at the recommended temperature, well closed and care is taken to prevent contamination during its use.

All other reagents are provided ready to use.

ADDITIONAL EQUIPMENT

- Moist chamber.
- Wash tray.
- Fluorescence microscope equipped with a 495 nm excitation filter and a 525 nm emission filter for FITC visualization.

SAMPLES

Serum or plasma collected by standard procedures. Stable 7 days at 2-8°C.

Dilute samples 1/20 in PBS (see Reagent Preparation) before the assay.

For titration of positive samples, make two-fold serial dilutions starting from 1/40 in PBS.

However, each laboratory should establish its own sample dilution based on the population characteristics.

PROCEDURE

1. Bring the reagents and samples to room temperature.
2. Place 1 drop (50 μ L) of the diluted sample or Control on each slide well, making sure that it is completely covered (Note 1).
3. Incubate the slide for 30 minutes at room temperature (15-30°C) into a moist chamber.
4. Drain sample drops off by gently tapping the inclined slide. Avoid cross-contamination of the sera.
5. Rinse gently the slide with PBS (see Reagent Preparation) (Note 2).
6. Wash thoroughly the slide by immersing in a washing tray filled with PBS for 5 minutes. Change PBS and repeat wash.
7. Carefully dry off the slides by using blotting paper. Keep the substrate moist along the procedure.

- Place 1 drop of Conjugate on each well. Incubate the slide for 30 minutes at room temperature (15-30°C) into a moist chamber.
- Wash (step 6) and dry (step 7).
- Place several drops of Mounting Medium on the slide and cover with a coverslip avoiding the formation of air bubbles.

READING

Examine the tissue section using the fluorescence microscope (250-400x). For best results, the slides should be read immediately, avoiding the center and the edge of the tissue.

Observation of specific fluorescent staining, described as follows at the recommended dilution, should be considered as a positive result.

ANA-homogenous: Homogeneous and uniform fluorescence throughout the nucleus of interphase cells. Strong fluorescence of mitotic cells.

AMA: Granular fluorescence of mitochondria in the cytoplasm of renal tubular cells.

ASMA: Staining of the muscularis mucosae, the muscle layers of the blood vessels and the interglandular fibers, in the rat stomach.

APCA: Reticular intracellular staining of parietal cells of rat gastric mucosa.

LKM: Type I. Strong staining of hepatocyte cytoplasm in liver and cytoplasm in proximal tubules in kidney. Negative distal tubules.

Positive sera may be tittered. The serum titer is defined as the highest dilution showing a positive result.

When none of the above specific stainings are observed, the result should be considered negative for these autoantibodies.

QUALITY CONTROL

Positive Control (CONTROL +) and Negative Control (CONTROL -) provided with kits cod 50065 should be tested together with the patients samples, in order to verify the assay performance.

Positive Control (CONTROL +) should give the above described specific staining.

Negative Control (CONTROL -) should not give any specific staining.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ASSAY CHARACTERISTICS

The IgG FITC/Evans conjugate is calibrated versus the International Standard of the WHO for sheep anti-human immunoglobulins conjugated with FITC.

The specificity of the ANA-Ho Positive Control has been verified against the AF/CDC1 reference serum from the Centers for Disease Control.

DIAGNOSTIC CHARACTERISTICS

ANA: Sensitivity of antinuclear antibodies determination is higher than 95% for systemic lupus erythematosus, although specificity is fairly low².

AMA: The presence of mitochondrial antibodies is associated with primary biliary cirrhosis in over than 95% of patients^{3,4}.

ASMA: Anti-smooth muscle antibodies are found in the sera of 52-85% of patients with autoimmune chronic active hepatitis and 22% of patients with primary biliary cirrhosis^{5,6}.

APCA: Antibodies to gastric parietal cells are found in 90% of patients with pernicious anemia, usually associated with other tissue specific autoimmune diseases⁷.

LKM: Anti-LKM type I antibodies are considered as markers of autoimmune hepatitis type II⁸.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

1. Avoid touching the tissue fixed into the wells along the procedure.
2. Use a squeeze bottle or a pipette to wash the slides, avoiding cross-contamination among the adjacent samples.

BIBLIOGRAPHY

1. Melnicoff MJ. Immunofluorescence Methods for Microscopic Analysis. In: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Hollingsworth PN et al. Antinuclear antibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
3. Leung PSC et al. Mitochondrial Autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
4. MacKay IR, and Gershwin ME. The autoantibodies of primary biliary cirrhosis: clinico pathological correlations. In: Van Venrooij WJ and Maini RN eds. Manual of Biological Markers of Diseases. Kluwer Academic Publishers, 1996.
5. Whittingham S and Mackay IR. Smooth muscle autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
6. Jacob G and Schoenfeld Y. Actin autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
7. Gleeson PA et al. Parietal Cell Autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
8. Homberg JC, et al. Chronic active hepatitis associated with antiliver/kidney microsome antibody type 1: a second type of "autoimmune" hepatitis. Hepatology 1987;7(6):1333-9.

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BIBLIOGRAPHIE

1. Melnicoff MJ. Immunofluorescence Methods for Microscopic Analysis. In: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Hollingsworth PN et al. Antinuclear antibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
3. Leung PSC et al. Mitochondrial Autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
4. MacKay IR, and Gershwin ME. The autoantibodies of primary biliary cirrhosis: clinico pathological correlations. In: Van Venrooij WJ and Maini RN eds. Manual of Biological Markers of Diseases. Kluwer Academic Publishers, 1996.
5. Whittingham S and Mackay IR. Smooth muscle autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
6. Jacob G and Schoenfeld Y. Actin autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
7. Gleeson PA et al. Parietal Cell Autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
8. Homberg JC, et al. Chronic active hepatitis associated with antiliver/kidney microsome antibody type 1: a second type of "autoimmune" hepatitis. *Hepatology* 1987;7(6):1333-9.

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BIBLIOGRAFIA

1. Melnickoff MJ. Immunofluorescence Methods for Microscopic Analysis. In: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Hollingsworth PN et al. Antinuclear antibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
3. Leung PSC et al. Mitochondrial Autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
4. MacKay IR, and Gershwin ME. The autoantibodies of primary biliary cirrhosis: clinico pathological correlations. In: Van Venrooij WJ and Maini RN eds. Manual of Biological Markers of Diseases. Kluwer Academic Publishers, 1996.
5. Whittingham S and Mackay IR. Smooth muscle autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
6. Jacob G and Schoenfeld Y. Actin autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
7. Gleeson PA et al. Parietal Cell Autoantibodies. In: James B. Peter and Yehuda Schoenfeld eds. Autoantibodies. Elsevier, 1996.
8. Homberg JC, et al. Chronic active hepatitis associated with antiliver/kidney microsome antibody type 1: a second type of "autoimmune" hepatitis. Hepatology 1987;7(6):1333-9.

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MANUFACTURER: BioSystems S.A.