ANTI-nDNA ANTIBODIES (nDNA)

Indirect Immunofluorescence CRITHIDIA LUCILIAE

A.MENARINI
diagnostics

COD 50063 20 x 12 tests

COD 50064 10 x 12 tests

Only for *in vitro* use in the clinical laboratory

INTENDED USE

Reagents for the qualitative determination of anti-nDNA antibodies in human serum by indirect immunofluorescence.

PRINCIPLE OF THE METHOD

Serum anti-nDNA antibodies bind to the corresponding antigen present in *Crithidia luciliae*. The resulting antigen-antibody complexes are detected by means of a fluorescein labeled anti-human immunoglobulin, and visualized with the aid of a fluorescence microscope¹.

CONTENTS

		COD 50063	COD 50064
SLIDE	Slides	10 x 6 tests	10 x 12 tests
PBS	PBS (10x)	1 x 100 mL	-
CONTROL +	nDNA Positive Control	1 x 1 mL	-
CONTROL -	Negative Control	1 x 1 mL	-
CONJUGATE	IgG FITC/Evans	1 x 5 mL	-
MOUNT	Mounting Medium	1 x 5 mL	-
COVER SLD	Coverslips	1 x 12	-

COMPOSITION

SLIDE *Crithidia luciliae* fixed on each well.

PBS Sodium phosphate 112.5 mmol/L, potassium phosphate 30 mmol/L, sodium chloride 1.15 mol/L, sodium azide 0.95 g/L, pH 7.2.

CONTROL+ Human serum containing anti-native DNA antibodies (nDNA), sodium azide 0.95 g/L.

CONTROL - Human serum, sodium azide 0.95 g/L.

CONJUGATE Goat anti-human IgG conjugated with fluorescein isothiocyanate (FITC), Evans blue 0.01 g/L, sodium azide 0.95 g/L.

MOUNT Mowiol 1%, Glycerol 87%, Tris 20 mmol/L, sodium azide 0.95 g/L.

Human sera used in the preparation of the positive and negative controls have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, the controls should be handled cautiously as potentially infectious.

STORAGE

Store at 2-8°C.

Reagents are stable until the expiry date shown on the label when stored tightly closed and if contamination is prevented during their use.

Indications of deterioration:

- Liquid components: Presence of particulate material, turbidity.
- Slides: rips in the sealing bag, macroscopic defects on the cell culture like scratches or monolayer peeling off.

AUXILIARY REAGENTS

PBS PBS (10x)

CONJUGATE IgG FITC/Evans

MOUNT Mounting Medium

CONTROL + nDNA Positive Control: Human serum containing anti-native DNA antibodies (nDNA), sodium azide 0.95 g/L.

CONTROL - Negative Control.

REAGENT PREPARATION

PBS: Dilute PBS 1/10 with distilled water. Stable 1 month at 2-8°C, if stored at the recommended temperature, well closed and care is taken to prevent contamination during its use.

All other reagents are provided ready to use.

ADDITIONAL EQUIPMENT

- Moist chamber
- Wash tray
- Fluorescence microscope equipped with a 495 nm excitation filter and a 525 nm emission filter for FITC visualization

SAMPLES

Serum or plasma collected by standard procedures. Stable 7 days at 2-8°C.

Dilute samples 1/10 in PBS (see Reagent Preparation) before the assay.

For titration of positive samples, make two-fold serial dilutions starting from 1/20 in PBS.


However, each laboratory should establish its own sample dilution based on the population characteristics.

PROCEDURE

1. Bring the reagents and samples to room temperature.
2. Place 1 drop (50 μ L) of the diluted sample or Control on each slide well, making sure that it is completely covered (Note 1).
3. Incubate the slide for 30 minutes at room temperature (15-30°C) into a moist chamber.
4. Drain sample drops off by gently tapping the inclined slide. Avoid cross-contamination of the sera.
5. Rinse gently the slide with PBS (see Reagent Preparation) (Note 2).
6. Wash thoroughly the slide by immersing in a washing tray filled with PBS for 5 minutes. Change PBS and repeat wash.
7. Carefully dry off the slides by using the blotting paper provided. Keep the substrate moist along the procedure.
8. Place 1 drop of Conjugate on each well. Incubate the slide for 30 minutes at room temperature (15-30°C) into a moist chamber.
9. Wash (step 6) and dry (step 7).
10. Place several drops of Mounting Medium on the slide and cover with a coverslip avoiding the formation of air bubbles.

READING

Examine the slide using the fluorescence microscope (250-400x). For best results, the slides should be read immediately. Select reading fields in the area indicated in the picture, between the center area and the edge area, with uniform spacing between cells. Fluorescent intensity in the center and edge areas is not representative of the slide preparation.

Crithidia luciliae is a hemoflagellate with a modified mitochondrion called kinetoplast containing double stranded DNA (nDNA), not associated to histones, which apparently lacks any other nuclear antigen. The kinetoplast is a rounded organelle smaller than the nucleus and located between this and the basal body (2,3,4).

Observation of specific fluorescent staining of the kinetoplasts at the recommended dilution should be considered as a positive result. Nuclei, basal body or flagellum might give fluorescent staining but only kinetoplast fluorescence must be scored.

Positive sera may be titered. The serum titer is defined as the highest dilution showing a positive result.

QUALITY CONTROL

Positive Control (CONTROL +) and Negative Control (CONTROL -) provided with kits cod 50063 should be tested together with the patients samples, in order to verify the assay performance.

Positive Control (CONTROL +) should give the above described specific staining.

Negative Control (CONTROL -) should not give any specific staining.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ASSAY CHARACTERISTICS

The IgG FITC/Evans conjugate is calibrated versus the International Standard of the WHO for sheep anti-human immunoglobulins conjugated with FITC.

DIAGNOSTIC CHARACTERISTICS

The *Crithidia luciliae* immunofluorescence test for anti-nDNA antibodies has a high diagnostic specificity, but a fairly high diagnostic sensitivity for Systemic Lupus Erythematosus (SLE). They are the most frequently detected autoantibodies associated to SLE: 95% in SLE patients with renal involvement, 50 – 70% in SLE patients without renal involvement and 40% in patients with inactive SLE. Anti-nDNA antibodies are rarely found in healthy individuals (5).

Menarini anti-nDNA antibodies kit was used to test 80 sera from SLE patients as well as healthy donors. The results are described as follows

	N	POS	NEG	Sens.	Spec.
SLE (<i>systemic lupus eritematosus</i>)	43	29	14	67%	100%
SLE without renal involvement	25	15	10	60%	100%
SLE with renal involvement	18	14	4	78%	100%
Other autoimmune diseases	20	0	20		
Healthy controls	17	0	17		

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

1. Avoid touching the cells fixed into the wells along the procedure.
2. Use a squeeze bottle or a pipette to wash the slides, avoiding cross-contamination among the adjacent samples.

BIBLIOGRAPHY

1. Melnicoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Aarden LA, de Groot ER and Feltkamp TEW. Immunology of DNA III. Crithidia luciliae, a simple substrate for the determination of anti-dsDNA with the immunofluorescence technique. Ann NY Acad Sci 1975; 254: 505-515.
3. Stingl G et al. An Immunofluorescence Procedure for the Demonstration of Antibodies to Native, Double-Stranded DNA and Circulating DNA-Anti-DNA Complexes. Clin Immunol and Immunopathol 1976; 6: 131-140.
4. Konstantinov KN and Russanova VR. Evidence for absence of histones in the Crithidia luciliae kinetoplast: a study with ant-H2A and monoclonal anti-H3 antibodies. British Journal of Dermatology 1987; 117: 451-456.
5. Conrad K, Schößler W, Hiepe F. Autoantibodies in Systemic Autoimmune Diseases. A Diagnostic Reference. Autoantigens, autoantibodies, autoimmunity. Vol. 2. Pabst Science Publishers, 2002.

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BIBLIOGRAFIA

1. Melnicoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Aarden LA, de Groot ER and Feltkamp TEW. Immunology of DNA III. Crithidia luciliae, a simple substrate for the determination of anti-dsDNA with the immunofluorescence technique. Ann NY Acad Sci 1975; 254: 505-515.
3. Stingl G et al. An Immunofluorescence Procedure for the Demonstration of Antibodies to Native, Double-Stranded DNA and Circulating DNA-Anti-DNA Complexes. Clin Immunol and Immunopathol 1976; 6: 131-140.
4. Konstantinov KN and Russanova VR. Evidence for absence of histones in the Crithidia luciliae kinetoplast: a study with ant-H2A and monoclonal anti-H3 antibodies. British Journal of Dermatology 1987; 117: 451-456.
5. Conrad K, Schöblier W, Hiepe F. Autoantibodies in Systemic Autoimmune Diseases. A Diagnostic Reference. Autoantigens, autoantibodies, autoimmunity. Vol. 2. Pabst Science Publishers, 2002.

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1. Melnickoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Aarden LA, de Groot ER and Feltkamp TEW. Immunology of DNA III. Crithidia luciliae, a simple substrate for the determination of anti-dsDNA with the immunofluorescence technique. Ann NY Acad Sci 1975; 254: 505-515.
3. Stingl G et al. An Immunofluorescence Procedure for the Demonstration of Antibodies to Native, Double-Stranded DNA and Circulating DNA-Anti-DNA Complexes. Clin Immunol and Immunopathol 1976; 6: 131-140.
4. Konstantinov KN and Russanova VR. Evidence for absence of histones in the Crithidia luciliae kinetoplast: a study with anti-H2A and monoclonal anti-H3 antibodies. British Journal of Dermatology 1987; 117: 451-456.
5. Conrad K, Schöblier W, Hiepe F. Autoantibodies in Systemic Autoimmune Diseases. A Diagnostic Reference. Autoantigens, autoantibodies, autoimmunity. Vol. 2. Pabst Science Publishers, 2002.

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BIBLIOGRAFIA

1. Melnicoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Aarden A, de Groot ER and Feltkamp TEW. Immunology of DNA III. Crithidia luciliae, a simple substrate for the determination of anti-dsDNA with the immunofluorescence technique. Ann NY Acad Sci 1975; 254: 505-515.
3. Stingl G et al. An Immunofluorescence Procedure for the Demonstration of Antibodies to Native, Double-Stranded DNA and Circulating DNA-Anti-DNA Complexes. Clin Immunol and Immunopathol 1976; 6: 131-140.
4. Konstantinov KN and Russanova VR. Evidence for absence of histones in the Crithidia luciliae kinetoplast: a study with anti-H2A and monoclonal anti-H3 antibodies. British Journal of Dermatology 1987; 117: 451-456.
5. Conrad K, Schöfler W, Hiepe F. Autoantibodies in Systemic Autoimmune Diseases. A Diagnostic Reference. Autoantigens, autoantibodies, autoimmunity. Vol. 2. Pabst Science Publishers, 2002.

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1. Melnickoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Aarden LA, de Groot ER and Feltkamp TEW. Immunology of DNA III. Crithidia luciliae, a simple substrate for the determination of anti-dsDNA with the immunofluorescence technique. Ann NY Acad Sci 1975; 254: 505-515.
3. Stingl G et al. An Immunofluorescence Procedure for the Demonstration of Antibodies to Native, Double-Stranded DNA and Circulating DNA-Anti-DNA Complexes. Clin Immunol and Immunopathol 1976; 6: 131-140.
4. Konstantinov KN and Russanova VR. Evidence for absence of histones in the Crithidia luciliae kinetoplast: a study with ant-H2A and monoclonal anti-H3 antibodies. British Journal of Dermatology 1987; 117: 451-456.
5. Conrad K, Schöblier W, Hiepe F. Autoantibodies in Systemic Autoimmune Diseases. A Diagnostic Reference. Autoantigens, autoantibodies, autoimmunity. Vol. 2. Pabst Science Publishers, 2002.

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LITERATUR

1. Melnickoff MJ. Immunofluorescence Methods for Microscopic Analysis. En: Howard GC ed. Methods in Nonradioactive Detection. Appleton & Lange, 1993.
2. Aarden LA, de Groot ER and Feltkamp TEW. Immunology of DNA III. Crithidia luciliae, a simple substrate for the determination of anti-dsDNA with the immunofluorescence technique. Ann NY Acad Sci 1975; 254: 505-515.
3. Stingl G et al. An Immunofluorescence Procedure for the Demonstration of Antibodies to Native, Double-Stranded DNA and Circulating DNA-Anti-DNA Complexes. Clin Immunol and Immunopathol 1976; 6: 131-140.
4. Konstantinov KN and Russanova VR. Evidence for absence of histones in the Crithidia luciliae kinetoplast: a study with ant-H2A and monoclonal anti-H3 antibodies. British Journal of Dermatology 1987; 117: 451-456.
5. Conrad K, Schöflner W, Hiepe F. Autoantibodies in Systemic Autoimmune Diseases. A Diagnostic Reference. Autoantigens, autoantibodies, autoimmunity. Vol. 2. Pabst Science Publishers, 2002.

MANUFACTURER: BioSystems S.A.
PRODUTTORE: Costa Brava 30
FABRICANTE: 08030 Barcelona (Spain)
FABRICANT: Tel: +34 93 311 00 00
ΚΑΤΑΣΚΕΥΑΣΤΗΣ: Fax: +34 93 311 76 09
HERSTELLER: Mail: Biosystems@biosystems.es

DISTRIBUTOR: A. Menarini Diagnostics SRL
DISTRIBUITORE: Via Sette Santi 3
DISTRIBUIDOR: 50131 Firenze
DISTRIBUTEUR: Tel. 0039-055-56801
ΔΙΑΝΟΜΕΑΣ: Fax 0039-055-5680-902
HÄNDLER: